


# Neurolyser XR Customer Onboarding Guide




## Site Survey: Pre-arrival



1

- **Overview** – Site provides FUSMobile information on imaging table and C-arm
- **Requirement:** C-arm video output type, side and bottom photos of table


## Installation: Day One



2

- **Overview:** : FUSMobile applications specialist arrives on site to install the Neurolyser XR and verify integration with the C-arm
- **Requirements:** Access to C-arm, rad-tech (if necessary), and cart
- **Duration:** 2 hours


## Clinical Training: Day One/Two



3

- **Overview:** Training presentation and hands on training is provided for physicians and supporting staff
- **Requirements:** Access to the C-arm and on-site staff (physicians, rad-techs, nurses, etc.)
- **Duration:** 30 minutes for presentation, 30 minutes for physicians, and 30 minutes for supporting staff

## Guided Procedures: Day Two/Three



4

- **Overview:** Clinicians schedule at least 3 procedures per doctor to be proctored by FUSMobile applications specialist
- **Requirements:** See patient selection guidelines
- **Duration:** 45 minutes to an hour per proctored patient
- **Note:** procedures can be done on either or both days